

# Barrel Shaker VKS 75 F control



## Description

The big heavy-duty shaker for barrels VKS 75 F control with touch display and programmable control has a horizontal shaking motion. With a shaking speed up to 80 rpm, a selected stroke between 120 - 150 mm, a loading capacity up to 75 kg and a respective rack system with 3 fixing elements the shaker is optimally suited for horizontal or vertical barrels up to 80 l capacity. **(floor installation only!)**

## Basic Equipment

Basic device **without shaking plate, without rack system**  
With touch display and programmable control

### **Loading Capacity**

Horizontal or vertical barrels up to 80 l capacity

[Loading details](#) against shaking speed.

### **Options:**

[USB 2.0 interface](#)

[RS 232 interface](#)

[Control interface](#)

[Safety Package](#)

## **Details**

<b>Motion:</b>	reciprocating (horizontal)
<b>Shaking platform:</b>	760 x 600 mm
<b>Max. load:</b>	75 kg
<b>Shaking speed:</b>	up to 80 rpm, in steps of 5
<b>Stroke:</b>	120 -150 mm in 10 mm steps
<b>Timer:</b>	programmable / continuous
<b>Order no.:</b>	6212 000
<b>Accessory:</b>	Shaking plate VKS (1-storey) Special rack system for barrels Special rack system for horizontal barrels

## **Technical specifications**

<b>Electrical supply:</b>	100 - 240 V, 50 / 60 Hz please indicate in case of order
<b>Enclosure protection:</b>	IP 21
<b>Heat emission:</b>	max. 15 W
<b>Ambient temperature:</b>	5°C to 50°C
<b>Relative humidity:</b>	~ 85 %
<b>Dimensions (w x d x h):</b>	1050 x 835 x 250 mm
<b>Weight:</b>	110 kg
<b>Total weight (Shaker + shaking plate + Rack system):</b>	132 kg