

Tube Roller RM 2



Description

Gentle tilting and rotating motion for uniform mixing

Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

Details

Motion:	rotating and tilting motion
Shaking platform:	200 x 300 mm
Max. load:	2 kg
Shaking speed:	5 - 80 rpm, in steps of 5
Angle of tilt:	3° fixed
Stroke at end of rolls:	15 mm

Dimensions of rolls: 300 x Ø 32 mm
Number of rolls: 6 pieces
Runtime: 0 - 99 h 59 min / continuous
Order No.: 6138 000
Accessory: Overhead Rack RM 2

Technical specifications

Electrical supply: 100 - 240 V, 50 / 60 Hz
Enclosure protection: IP 21
Heat emission: max. 5 W
Ambient temperature: 5°C to 45°C
Relative humidity: ~ 85 % (non-condensing)
Dimensions (w x d x h): 244 x 465 x 115 mm
Weight: 7,3 kg
